



Test Catalog

Diagnostic. Prognostic. Predictive. Predisposition.



NeoTYPE® Thyroid Profile

Alternative Name

Thyroid Profile

Methodology

Molecular

Test Description

The NeoTYPE Thyroid Tumor Profile analyzes 32 biomarkers through a combination of next-generation sequencing (NGS), FISH, and IHC as listed below. Test orders include summary interpretation of all results to help guide treatment decisions. If Pan-TRK IHC is expressed or equivocal, reflex to either NTRK NGS Fusion Panel (Default) or NTRK 1-3 FISH Panel will be added. A microsatellite instability (MSI) NGS result of “indeterminate” will create a reflex to MSI by PCR as long as the tumor percentage is ≥40% and paired normal tissue is available.

- NGS (26 genes + 2 biomarkers): AKT1, ALK, ARID1A, ATM, BRAF, CDKN2A, CTNNB1, ERBB2, ERBB4, HRAS, KRAS, MEN1, MET, Microsatellite Instability (MSI), NF1, NF2, NRAS, PIK3CA, PTEN, RET, SMAD4, SMO, SRC, TERT Promoter, TP53, TSC1, TSC2, Tumor Mutation Burden (TMB)
- FISH (2 FISH): MET, RET (tech-only available)
- IHC (2 biomarkers): PD-L1 LDT, Pan-TRK (tech-only available for PD-L1)

Clinical Significance

The NeoTYPE Thyroid Profile is useful to classify fine needle aspirates of thyroid nodules that are indeterminate or suspicious on cytology. Presence of mutations or gene rearrangements as detected by FISH predicts malignancy with varying degrees of specificity depending on the gene mutated and histological subtype. BRAF mutation V600E is associated with poor prognosis in papillary thyroid carcinoma (PTC).

Specimen Requirements

- **FFPE tissue:** Paraffin block preferred. Please use 10% buffered formalin fixative. Do not use zinc fixatives.

Storage & Transportation

Use cold pack for transport, making sure cold pack is not in direct contact with specimen. All slides can be packed at room temperature.

CPT Code(s)*

81445x1, 88377x2, 88360x1, 88342x1; add 81479x1 if reflexed to NTRK NGS Fusion Panel (default) or 88374x3 automated (88377x3 manual) if reflexed to NTRK 1-3 FISH Panel

Medicare MoIDX CPT Code(s)*

81479

New York Approved

Yes

Level of Service

Global

Turnaround Time

14 days

NeoGenomics Laboratories is a specialized oncology reference laboratory providing the latest technologies, testing partnership opportunities, and interactive education to the oncology and pathology communities. We offer the complete spectrum of diagnostic services in molecular testing, FISH, cytogenetics, flow cytometry, and immunohistochemistry through our nation-wide network of CAP-accredited, CLIA-certified laboratories.

Committed to research as the means to improve patient care, we provide Pharma Services for pharmaceutical companies, in vitro diagnostic manufacturers, and academic scientist-clinicians. We promote joint publications with our client physicians. NeoGenomics welcomes your inquiries for collaborations. Please contact us for more information.

*The CPT codes provided with our test descriptions are based on AMA guidelines and are for informational purposes only. Correct CPT coding is the sole responsibility of the billing party.

Please direct any questions regarding coding to the payor being billed.



9490 NeoGenomics Way
Fort Myers, FL 33912
Phone: 239.768.0600/ Fax: 239.690.4237
neogenomics.com

© 2024 NeoGenomics Laboratories, Inc. All Rights Reserved.
All other trademarks are the property of their respective owners
Rev. 052024